# Highlights of the Z5 Bone Level implant

The world's first 100% ceramic, two-piece, screw-retained implant-abutment connection

### **Z5-BL 40-10**

> Implant 4 mm × 10 mm

#### Aesthetics for professional care

- Wide range of straight and angled abutments for optimal esthetic results and flexible functionality
- A wide range of options for designing the exit profile

#### 100% ceramic

- 100 % ceramic, screw-retained, two-piece bone level implant
- Surgical procedure analogous to two-piece titanium implantology
- Closed healing, augmentation and simultaneous sinus lift is possible in suitable cases

## Experience from over

60,000 implants placed

- Established ennosal design, based on tissue-level predecessors, ensures transferable results
- Self-tapping thread with bone chip reservoir reduces insertion torque
- Excellent osseointegration through Z-Systems's proven, patented proprietary SLM® laser surface technology and gentle plasma sterilization

#### A Unique, strong connection

- The conical screw connection with internal thread for microgap-free connection and prevention of a pump effect.
- **NEW:** ceramic abutment screw for 100% all-ceramic restoration
- A titanium screw is also available as an option
- Platform switch for a more stable and enlarged mucosal cuff and for reduced loading of the crestal bone

